

ERCC6L Rabbit Polyclonal Antibody



CAB14137

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

141kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54821

Uniprot

Q2NKX8

Synonyms:

ERCC6L; PICH; RAD26L

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ERCC6L (NP_060139.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
