
Product Information

Product SKU:	CAB14139	Gene ID:	55072	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human, Mouse

Additional Information

Observed MW:	120kDa	Conjugate:	Unconjugated
Calculated MW:	120kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. The encoded protein is the E3 ubiquitin-protein ligase component of the linear ubiquitin chain assembly complex. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:100

Synonyms: HOIP; Paul; ZIBRA; RNF31

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human RNF31 (NP_060469.4).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.