

# PLCD4 Rabbit Polyclonal Antibody



CAB14154

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

31kDa/87kDa

### Applications:

WB

### Reactivity:

Mouse

## Protein Background

This gene encodes a member of the delta class of phospholipase C enzymes. Phospholipase C enzymes play a critical role in many cellular processes by hydrolyzing phosphatidylinositol 4, 5-bisphosphate into two intracellular second messengers, inositol 1, 4, 5-trisphosphate and diacylglycerol. Expression of this gene may be a marker for cancer.

## Immunogen information

### Gene ID:

84812

### Uniprot

Q9BRC7

### Synonyms:

PLCD4

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PLCD4 (NP\_116115.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---