

# GSTO2 Rabbit Polyclonal Antibody



CAB14156

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

24kDa/25kDa/28kDa

### Applications:

IF

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

IF 1:50 - 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

The protein encoded by this gene is an omega class glutathione S-transferase (GST). GSTs are involved in the metabolism of xenobiotics and carcinogens. Four transcript variants encoding different isoforms have been found for this gene.

## Immunogen information

### Gene ID:

119391

### Uniprot

Q9H4Y5

### Synonyms:

GSTO2; GSTO 2-2; bA127L20.1; GSTO2-2

### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 124-243 of human GSTO2 (NP\_899062.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunofluorescence analysis of U2OS cells using GSTO2 antibody (CAB14156). Blue: DAPI for nuclear staining.