

PLCXD2 Rabbit Polyclonal Antibody



CAB14164

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

34kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

257068

Uniprot

Q0VAA5

Synonyms:

PLCXD2


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human PLCXD2 (NP_001172035.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of A549 cells using PLCXD2 antibody (CAB14164).