

Product Information

Product SKU:	CAB14167	Gene ID:	339967	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable serine-type endopeptidase activity. Predicted to be involved in proteolysis. Predicted to be located in extracellular region. Predicted to be integral component of plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	HESP; ECRG1; HATL1; TMPRSS11A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 300-421 of human TMPRSS11A (NP_872412.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.