

CAB14197

Product Information

Product SKU:	CAB14197	Gene ID:	128239	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	170kDa	Conjugate:	Unconjugated
Calculated MW:	185kDa	Isotype:	IgG

Immunogen Information

Background: Enables calmodulin binding activity and myosin VI light chain binding activity. Predicted to be involved in regulation of actin cytoskeleton organization. Predicted to act upstream of or within several processes, including intracellular signal transduction; positive regulation of macromolecule metabolic process; and positive regulation of mammary gland epithelial cell proliferation. Predicted to be located in cytosol. Predicted to be active in cytoplasm.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: IQGAP3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1437-1631 of human IQGAP3 (NP_839943.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.