

OPRD1 Rabbit Polyclonal Antibody



CAB14207

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

40kDa

Calculated MW:

40kDa

Applications:

WB

Reactivity:

Human, Mouse

Protein Background

Immunogen information

Gene ID:

4985

Uniprot

P41143

Synonyms:

OPRD1; OPRD

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 283-372 of human OPRD1 (NP_000902.3).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using OPRD1 antibody (CAB14207) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.