

WDR33 Rabbit Polyclonal Antibody
CAB14208



Product Information

Product SKU:	CAB14208	Gene ID:	74320	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	145kDa	Conjugate:	Unconjugated
Calculated MW:	30kDa/38kDa/145kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in mRNA polyadenylation. Predicted to act upstream of or within mRNA processing. Located in nucleus. Orthologous to human WDR33 (WD repeat domain 33).

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	WDC146; 1110001N06Rik; 2310011G05Rik; 2810021O11Rik; 8430413N20Rik; WDR33
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of mouse WDR33 (NP_001164437.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.