

CAB14209

Product Information

Product SKU:	CAB14209	Gene ID:	50809	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables DNA binding activity and nucleosome binding activity. Involved in several processes, including cellular response to hypoxia; heterochromatin organization; and regulation of nucleus size. Located in chromosome and nuclear speck.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	HP1BP74; HP1-BP74; HP1BP3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human HP1BP3 (NP_057371.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.