

CAB14215

Product Information

Product SKU:	CAB14215	Gene ID:	54097	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	26kDa	Conjugate:	Unconjugated
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

Background: Involved in insulin secretion. Located in extracellular exosome.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	2-21; ORF9; PANDER; PRED44; C21orf11; C21orf76; FAM3B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human FAM3B (NP_478066.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.