

CAB14218

Product Information

Product SKU:	CAB14218	Gene ID:	24139	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	71kDa	Conjugate:	Unconjugated
Calculated MW:	71kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables microtubule binding activity. Involved in negative regulation of microtubule polymerization and regulation of microtubule nucleation. Colocalizes with mitotic spindle.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ELP70; EMAP2; EMAP-2; EML2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 150-415 of human EML2 (NP_036287.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.