

CCDC124 Rabbit Polyclonal Antibody
CAB14298



Product Information

Product SKU:	CAB14298	Gene ID:	115098	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

Background:	Enables RNA binding activity. Predicted to be involved in cell division. Located in cytosol and plasma membrane.
Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:100
Synonyms:	Lso2; oxs1; CCDC124
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human CCDC124 (NP_612451.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.