
Product Information

Product SKU:	CAB14303	Gene ID:	254170	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	63kDa	Conjugate:	Unconjugated
Calculated MW:	63kDa	Isotype:	IgG

Immunogen Information

Background: This locus represents an member of the F-box gene family. The encoded protein contains an F-box motif and a domain that might form a structure similar to a leucine-rich repeat found in placental RNase inhibitor. This locus may be associated with copy number variation of UGT2B17 (GeneID 7367), which has been associated with susceptibility to osteoporosis.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: Fbx33; BMND12; c14_5247; FBXO33

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 346-555 of human FBXO33 (NP_976046.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.