

Product Information

Product SKU:	CAB14310	Gene ID:	51504	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein heterodimerization activity and protein methyltransferase activity. Involved in macromolecule methylation and positive regulation of rRNA processing. Located in nucleoplasm and perinuclear region of cytoplasm. Part of protein-containing complex.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: TRM112; HSPC152; HSPC170; hTrm112; TRMT11-2; TRMT112

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human TRMT112 (NP_057488.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.