

RFTN2 Rabbit Polyclonal Antibody



CAB14340

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

55kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

130132

Uniprot

Q52LD8

Synonyms:

RFTN2; C2orf11; Raftlin-2; raftlin-2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 412-501 of human RFTN2 (NP_653230.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
