

Product Information

Product SKU:	CAB14353	Gene ID:	54532	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	120kDa	Conjugate:	Unconjugated
Calculated MW:	121kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to enable thiol-dependent deubiquitinase. Predicted to be involved in response to auditory stimulus and sensory perception of sound. Predicted to act upstream of or within action potential and neuron apoptotic process. Predicted to be located in bicellular tight junction. Predicted to be active in cell-cell junction.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PFIC7; USP53
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 886-1073 of human USP53 (NP_061923.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.