

Product Information

Product SKU:	CAB14383	Gene ID:	7579	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	117kDa	Conjugate:	Unconjugated
Calculated MW:	118kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be located in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	KOX29; ZNF31; ZFP-31; ZNF360; ZSCAN20
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ZSCAN20 (NP_660281.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.