

# RPUSD2 Rabbit Polyclonal Antibody



CAB14397

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

54kDa/61kDa

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

27079

#### Uniprot

Q8IZ73

#### Synonyms:

RPUSD2; C15orf19; C18B11

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 286-545 of human RPUSD2 (NP\_689473.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---