

TEKT2 Rabbit Polyclonal Antibody



CAB14398

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

49kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene product belongs to the tektin family of proteins. Tektins comprise a family of filament-forming proteins that are coassembled with tubulins to form ciliary and flagellar microtubules. This gene is expressed in the testis and its protein is localized to the flagella of the sperms, indicating that it may play a role in spermatogenesis.

Immunogen information

Gene ID:

27285

Uniprot

Q9UIF3

Synonyms:

TEKT2; TEKTB1; TEK TIN-T; h-tektin-t; tektin-2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-430 of human TEKT2 (NP_055281.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
