

ZNF585A Rabbit Polyclonal Antibody



CAB14457

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

46kDa/87kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

199704

Uniprot

Q6P3V2

Synonyms:

ZNF585A

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human ZNF585A (NP_001275729.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
