

ZFP82 Rabbit Polyclonal Antibody



CAB14461

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

62kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

284406

Uniprot

Q8N141

Synonyms:

ZFP82; ZNF545

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 50-170 of human ZFP82 (NP_597723.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
