

RPUSD3 Rabbit Polyclonal Antibody



CAB14462

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

14kDa/36kDa/38kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

285367

Uniprot

Q6P087

Synonyms:

RPUSD3

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 32-351 of human RPUSD3 (NP_775930.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
