

C7orf50 Rabbit Polyclonal Antibody



CAB14484

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

22kDa

Applications:

WB IHC

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

84310

Uniprot

Q9BRJ6

Synonyms:

C7orf50; YCR016W

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human C7orf50 (NP_115726.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
