

CAB14492

Product Information

Product SKU:	CAB14492	Gene ID:	6623	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	16kDa	Conjugate:	Unconjugated
Calculated MW:	13kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the synuclein family of proteins which are believed to be involved in the pathogenesis of neurodegenerative diseases. Mutations in this gene have also been associated with breast tumor development.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	SR; BCSG1; gamma-Synuclein
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 50-127 of human gamma-Synuclein (NP_003078.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.