

CAB14518

Product Information

Product SKU:	CAB14518	Gene ID:	84892	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	67kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: AGO61; GTDC2; MDDGA8; MDDGC8; C3orf39; POMGNT2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 160-340 of human POMGNT2 (NP_116195.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.