

POMGNT2 Rabbit Polyclonal Antibody



CAB14518

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

75kDa

Calculated MW:

66kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

Immunogen information

Gene ID:

84892

Uniprot

Q8NAT1

Synonyms:

POMGNT2; AGO61; C3orf39; GTDC2; MDDGA8


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 160-340 of human POMGNT2 (NP_116195.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using POMGNT2 antibody (CAB14518) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.