

ZDHHC14 Rabbit Polyclonal Antibody



CAB14538

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

51kDa/53kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79683

Uniprot

Q8IZN3

Synonyms:

ZDHHC14; NEW1CP

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 273-488 of human ZDHHC14 (NP_078906.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
