

# NCR1 Rabbit Polyclonal Antibody



CAB14544

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

21kDa/22kDa/23kDa/32kDa/  
34kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

9437

#### Uniprot

O76036

#### Synonyms:

CD335; LY94; NK-p46; NKP46; NCR1

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 151-303 of human NCR1 (NP\_001138929.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---