

# FAM13A Rabbit Polyclonal Antibody



CAB14576

---

## Product Information

<b>Product SKU:</b>	CAB14576	<b>Gene ID:</b>	10144	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

## Additional Information

<b>Observed MW:</b>	116kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	117kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Predicted to enable GTPase activator activity. Predicted to be involved in regulation of small GTPase mediated signal transduction. Predicted to be located in cytosol. Implicated in chronic obstructive pulmonary disease.

<b>Recommended Dilution:</b>	WB,1:500 - 1:2000
<b>Synonyms:</b>	FAM13A1; ARHGAP48; FAM13A
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human FAM13A (NP_055698.2).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.