

**CAB14589**

## Product Information

<b>Product SKU:</b>	CAB14589	<b>Gene ID:</b>	54850	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

## Additional Information

<b>Observed MW:</b>	37kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	37kDa	<b>Isotype:</b>	IgG

## Immunogen Information

**Background:** Members of the F-box protein family, such as FBXL12, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** Fbl12; FBXL12

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human FBXL12 (NP\_060173.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.