

Product Information

Product SKU:	CAB14595	Gene ID:	64123	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	78kDa	Conjugate:	Unconjugated
Calculated MW:	78kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable G protein-coupled receptor activity. Predicted to be involved in adenylate cyclase-activating G protein-coupled receptor signaling pathway. Predicted to be located in cytoplasmic vesicle and plasma membrane. Predicted to be integral component of plasma membrane. Biomarker of glioblastoma and hypertrophic cardiomyopathy.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ETL; ELTD1; KPG_003; ADGRL4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 110-406 of human ADGRL4 (NP_071442.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.