

DHRS11 Rabbit Polyclonal Antibody



CAB14597

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

20kDa/28kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79154

Uniprot

Q6UWP2

Synonyms:

DHRS11; ARPG836; SDR24C1; spDHRS11

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DHRS11 (NP_077284.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
