

SLR1 Rabbit Polyclonal Antibody



CAB14614

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

65.4kDa

Calculated MW:

65kDa

Applications:

WB

Reactivity:

Oryza sativa

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4333860

Uniprot

Synonyms:

SLR1


Immunogen:

Recombinant protein of Oryza sativa SLR1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - SLR1 Rabbit pAb (CAB14614)