

CAB14646

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB14646 | Gene ID: | 3150 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Mouse,Rat |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | - | Conjugate: | Unconjugated |
| Calculated MW: | 11kDa | Isotype: | IgG |

Immunogen Information

Background: The protein encoded by this gene binds nucleosomal DNA and is associated with transcriptionally active chromatin. Along with a similar protein, HMG17, the encoded protein may help maintain an open chromatin configuration around transcribable genes.

| | |
|-----------------------|---|
| Recommended Dilution: | IHC-P,1:50 - 1:200 |
| Synonyms: | HMG14; HMGN1 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human HMGN1 (NP_004956.5). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |