

IFT57 Rabbit Polyclonal Antibody



CAB14652

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

49kDa

Applications:

WB IHC

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55081

Uniprot

Q9NWB7

Synonyms:

IFT57; ESRRBL1; HIPPI; MHS4R2


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 170-429 of human IFT57 (NP_060480.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human uterine cancer using IFT57 antibody (CAB14652) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat ovary using IFT57 antibody (CAB14652) at dilution of 1:100 (40x lens).