

CAB14661

Product Information

Product SKU:	CAB14661	Gene ID:	54855	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA adenylyltransferase activity. Involved in mRNA stabilization. Located in cytoplasm and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FAM46C
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human FAM46C (NP_060179.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.