

ZNF83 Rabbit Polyclonal Antibody



CAB14674

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

59kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55769

Uniprot

P51522

Synonyms:

ZNF83; HPF1; ZNF816B

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZNF83 (NP_060770.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
