

Product Information

Product SKU:	CAB14686	Gene ID:	55852	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	130kDa	Conjugate:	Unconjugated
Calculated MW:	125kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable lipid binding activity. Predicted to be involved in lipid transport. Predicted to be located in endoplasmic reticulum membrane and nuclear membrane. Predicted to be integral component of membrane. Predicted to be active in endoplasmic reticulum.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	HT008; TMEM96; TEX2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 915-1134 of human TEX2 (NP_060939.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.