

Product Information

Product SKU:	CAB14733	Gene ID:	2637	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	37kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in regulation of nervous system development and regulation of transcription by RNA polymerase II. Predicted to act upstream of or within several processes, including branching involved in blood vessel morphogenesis; nervous system development; and neural crest cell migration. Predicted to be part of chromatin. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	GBX2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human GBX2 (NP_001476.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.