

Product Information

Product SKU:	CAB14759	Gene ID:	5067	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	113kDa	Conjugate:	Unconjugated
Calculated MW:	113kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in cell adhesion and nervous system development. Predicted to be located in extracellular region and plasma membrane. Predicted to be active in neuron projection.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PCS; PANG; BIG-1; CNTN3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human CNTN3 (NP_065923.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.