

Product Information

Product SKU:	CAB14767	Gene ID:	5365	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	172kDa	Conjugate:	Unconjugated
Calculated MW:	207kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the plexin family. It functions as a receptor for semaphorin 5A, and plays a role in axon guidance, invasive growth and cell migration. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PLEXR; PLXN6; PLEXB3; PLXNB3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1050-1250 of human PLXNB3 (NP_001156729.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.