

CAB14792

Product Information

Product SKU:	CAB14792	Gene ID:	7294	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	85kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable non-membrane spanning protein tyrosine kinase activity. Involved in positive regulation of interferon-gamma-mediated signaling pathway; positive regulation of macromolecule metabolic process; and protein autophosphorylation. Located in cytoplasm and nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	RLK; TKL; BTKL; PTK4; PSCTK5; TXK
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TXK (NP_003319.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.