

---

## Product Information

<b>Product SKU:</b>	CAB14803	<b>Gene ID:</b>	7730	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

## Additional Information

<b>Observed MW:</b>	55kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	55kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in negative regulation of transcription by RNA polymerase II. Located in blood microparticle.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** PIGX; ZNF177

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 50-110 of human ZNF177 (NP\_003442.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.