

ZNF200 Rabbit Polyclonal Antibody



CAB14804

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

45kDa

Applications:

IHC

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7752

Uniprot

P98182

Synonyms:

ZNF200

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-395 of human ZNF200 (NP_003445.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
