

SLC5A6 Rabbit Polyclonal Antibody



CAB14814

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

68kDa

Applications:

IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

8884

Uniprot

Q9Y289

Synonyms:

SLC5A6; SMVT


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human SLC5A6 (NP_066918.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human liver using SLC5A6 antibody (CAB14814) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat heart using SLC5A6 antibody (CAB14814) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse heart using SLC5A6 antibody (CAB14814) at dilution of 1:100 (40x lens).