

KIAA0513 Rabbit Polyclonal Antibody
CAB14829



Product Information

Product SKU:	CAB14829	Gene ID:	9764	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	47kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be located in cytoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	KIAA0513
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human KIAA0513 (NP_001273494.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.