

ZNF273 Rabbit Polyclonal Antibody



CAB14846

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

58kDa/64kDa

Applications:

IHC

Reactivity:

Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene is a member of the krueppel C2H2-type zinc-finger protein family and encodes a protein with 13 C2H2-type zinc fingers and a KRAB domain. This nuclear protein is involved in transcriptional regulation. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

10793

Uniprot

Q14593

Synonyms:

ZNF273; HZF9

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 210-260 of human ZNF273 (NP_066971.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
