
Product Information

Product SKU:	CAB14858	Gene ID:	23177	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	81kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein domain specific binding activity and protein kinase binding activity. Involved in centriole-centriole cohesion and protein localization to organelle. Located in several cellular components, including centriolar satellite; cytosol; and nucleoplasm.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: KIAA0582; CEP68

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human CEP68 (NP_001306030.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.