

# LRP10 Rabbit Polyclonal Antibody



CAB14867

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

60kDa/76kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

26020

#### Uniprot

Q7Z4F1

#### Synonyms:

LRP10; LRP-10; LRP9; MST087; MSTP087

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 17-120 of human LRP10 (NP\_054764.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---